

C-sharp windows example

Online version: https://wiki.advacam.cz/wiki/C-sharp_windows_example



Contents

Overview & notes	3
Window	3
Code	4
Related	10

Overview & notes

This is C# example of the Windows program with list devices, simple measurement and view of the output data.

Timepix3 only.


Notes:

- Use the release/64 bit configuration
- The working directory is directory with the exe file. Typically **project\bin\Release**. Copy pixet.ini and other auxilliary files here.
- In the MS Visual studio 2019, **project first not working**. You must click Properties, change .NET version to old, save it, change .NET version back to actual and save. Now project can work.

Window

Create window with components:

1. lstMsg - Listbox for messages
2. lstDevs - Listbox for devices list
3. btnInit - Button for init the Pixet core
4. btnClearMsg - Button for clear lstMsg
5. btnExit - Button for exit the Pixet Core
6. pbFrame1 - PictureBox for first output subframe, size equal 256x256 inside, Border Fixed3D
7. pbFrame2 - PictureBox for second output subframe, size equal 256x256 inside, Border Fixed3D
8. txtMeasTime - Textbox for measure time input, default text "1"
9. rbOPM0 - Radiobutton for operation mode 0
10. rbOPM2 - Radiobutton for operation mode 2
11. checkLogView - Checkbox to enable logaritimized view
12. checkColView - Checkbox to enable colorized view
13. btnMeasSinTpx3 - Button to run the pxcMeasureSingleFrameTpx3 example measuring
14. labImgInfo - Label for informations about measured images



C-sharp windows example - window edit screenshot

Code

```
using System;
using System.Drawing;
using System.Linq;
using System.Text;
using System.Windows.Forms;
using System.Runtime.InteropServices;

namespace WindowsFormsApp1 {
    public partial class Form1 : Form {

        [DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
        public static extern int pxcInitialize(Int32 a, UInt64 b);
    }
}
```

```

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcExit();

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcGetDevicesCount();

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcGetDeviceChipID(UInt32 deviceIndex, UInt32 chipIndex,
StringBuilder chipIDBuffer, UInt32 size);

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcGetDeviceName(UInt32 deviceIndex, StringBuilder
nameBuffer, UInt32 size);

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcGetDeviceChipCount(UInt32 deviceIndex);

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcGetDeviceInfo(UInt32 deviceIndex, UInt64 devInfo); //
devInfo is struct, but not used here

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcMeasureSingleFrameTpx3(UInt32 deviceIndex, double
frameTime, [Out] double[] frameToaITot, [Out] UInt16[] frameTotEvent, ref UInt32 size,
UInt32 trgStg = 0);
[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcGetDeviceDimensions(UInt32 deviceIndex, ref UInt32
width, ref UInt32 height);

[DllImport("pxcore.dll", CallingConvention = CallingConvention.Cdecl)]
public static extern int pxcSetTimepix3Mode(UInt32 deviceIndex, int mode);

int devicesCout = 0;
int[] devTypesN = Enumerable.Repeat(-1, 1000).ToArray();
UInt32 deviceIdx = uint.MaxValue;
UInt32 devXsize = 256, devYsize =256;
UInt32 devPixels = 65536;

public Form1() {
    InitializeComponent();
}

void msgToList(String msg) { //
=====
    lstMsg.Items.Add(msg);
}
void msgToList(String msg, int rc) {
    lstMsg.Items.Add($"{msg:S} rc={rc:D} (0 is OK)");
}

```

```

void msgToList(String msg1, int n, String msg2) {
    lstMsg.Items.Add($"{msg1:S} {n:D} {msg2:S}");
}
void msgToList(String msg1, UInt32 n, String msg2) {
    lstMsg.Items.Add($"{msg1:S} {n:D} {msg2:S}");
}

int setPixelRGB(int r, int g, int b) { //
=====
    return b + (g << 8) + (r << 16);
}

// value 0-1280 to RGB color or 0-255 to grayscale
int valToRGB(double val) { //
=====
    int v = (int)val;

    if (!checkColView.Checked) return setPixelRGB(v, v, v);

    if (v < 256) return setPixelRGB(0, 0, v); // 0 black to blue
    v -= 256;
    if (v < 256) return setPixelRGB(0, v, 255); // 256 blue to blue-green
    v -= 256;
    if (v < 256) return setPixelRGB(0, 255, 255 - v); // 512 blue-green to green
    v -= 256;
    if (v < 256) return setPixelRGB(v, 255, 0); // 768 green to yellow
    v -= 256;
    if (v < 256) return setPixelRGB(255, 255 - v, 0); // 1024 yellow to red
    v -= 256;
    if (v < 256) return setPixelRGB(255, 0, v); // 1280 red to violet
    v -= 256;
    if (v < 256) return setPixelRGB(255, v, 255); // 1536 violet to white
    else return setPixelRGB(255, 255, 255); // over 1792
}

const double cColorizeIniMin = 1e30, cColorizeIniMax = -1e30;
double colorizeMin = cColorizeIniMin, colorizeMax = cColorizeIniMax;
int colorizePixCount = 0;
// normalize data (, logarithmize), convert to color RGB (or grayscale), count non-
zero, count values
// note: logarithmize not usable if data<=0
void colorizeData(double[] dataIn, int[] dataOut) { //
=====
    double colMin=0, colMax=0, pixVal, p;
    colorizeMin = cColorizeIniMin;
    colorizeMax = cColorizeIniMax;
    colorizePixCount = 0;

    for (int n = 0, i = 0; n < devPixels; n++) { // find min, max

```

```

        if (dataIn[n] != 0) {
            if (checkLogView.Checked) pixVal = Math.Log(dataIn[n]); else pixVal =
dataIn[n];
            if (dataIn[n] > colorizeMax) { colorizeMax = dataIn[n]; colMax =
pixVal; };
            if (dataIn[n] < colorizeMin) { colorizeMin = dataIn[n]; colMin =
pixVal; };
            colorizePixCount++;
        }
    }

    if (checkColView.Checked) p = 1792.0 / (colMax - colMin);
    else p = 255.0 / (colMax - colMin);
    if (p < 0) p = -p;

    for (int n = 0; n < devPixels; n++) { // colorize data
        if (dataIn[n] != 0) {
            if (checkLogView.Checked) pixVal = Math.Log(dataIn[n]); else pixVal =
dataIn[n];
            dataOut[n] = valToRGB((pixVal - colMin) * p);
        } else {
            dataOut[n] = setPixelRGB(128, 128, 128);
        }
    }
}

void colorizeData(UInt16[] dataIn, int[] dataOut) { //
=====
    double[] tmpColData = new double[devPixels];
    for (int n = 0; n < devPixels; n++) tmpColData[n] = (double)dataIn[n];
    colorizeData(tmpColData, dataOut);
}

private void btnInit_Click(object sender, EventArgs e) { //
=====
    btnInit.Enabled = false;
    int rc = pxcInitialize(0, 0);
    msgToList("pxcInitialize", rc);
    rc = pxcGetDevicesCount();
    msgToList("pxcGetDevicesCount", rc, "(>=0 devCnt, <0 error)");
    lstDevs.Items.Clear();
    if (rc>0) {
        devicesCout = rc;
        for (UInt32 n=0; n<devicesCout; n++) {
            UInt32 bufLen = 100;
            StringBuilder buff = new StringBuilder((int)bufLen);

            rc = pxcGetDeviceName(n, buff, bufLen);
            String devName;
            if (rc == 0) { devName = buff.ToString(); }

```

```

else { devName = "(failed)"; }

rc = pxcGetDeviceChipID(n, 0, buff, bufLen);
String chipID;
if (rc == 0) { chipID = buff.ToString(); }
else { chipID = "(failed)"; }

int chipCnt = pxcGetDeviceChipCount(n);

String[] devTypes = { "unknown", "TPX", "MPX3", "TPX3", "TPX2" };
rc = pxcGetDeviceInfo(n, 0);
String dt;
if (rc >= 0 && rc < devTypes.Length) {
    dt = devTypes[rc];
    devTypesN[n] = rc;
} else { dt = "(failed)"; }

lstDevs.Items.Add($"{n:D}: ID:{devName:S}, ID0:{chipID:S},
chipCnt:{chipCnt:D}, type:{dt:S}");
}
btnExit.Enabled = true;
} else {
    msgToList("No devs found");
    btnInit.Enabled = true;
    btnExit.Enabled = false;
    btnMeasSinTpx3.Enabled = false;
}
}

private void btnClearMsg_Click(object sender, EventArgs e) { //
=====
    lstMsg.Items.Clear();
}

private void Form1_Load(object sender, EventArgs e) {
    lstMsg.Items.Add("Messages");
    lstDevs.Items.Add("Devices list - click the 'Init core' button");
}

private void Form1_FormClosing(object sender, FormClosingEventArgs e) {
    pxcExit();
}

private void btnExit_Click(object sender, EventArgs e) { //
=====
    int rc = pxcExit();
    msgToList("pxcExit", rc);
    lstDevs.Items.Clear();
    deviceIdx = uint.MaxValue;
}

```



```

    devicesCout = 0;

    btnInit.Enabled = true;
    btnExit.Enabled = false;
    btnMeasSinTpx3.Enabled = false;
}

private void lstDevs_SelectedIndexChanged(object sender, EventArgs e) { //
=====
    UInt32 i = (UInt32)lstDevs.SelectedIndex;
    if (devTypesN[i]==3) {
        deviceIdx = i;
        msgToList("Selected device:", deviceIdx, "");
        UInt32 width = 256, height = 256;
        int rc = pxcGetDeviceDimensions(deviceIdx, ref width, ref height);
        if (rc==0) {
            devXsize = width; devYsize = height;
            devPixels = devXsize * devYsize;
            msgToList($"Dev dimm: w:{devXsize:D}, h:{devYsize:D}");
        } else {
            msgToList("pxcGetDeviceDimensions", rc);
        }
        btnMeasSinTpx3.Enabled = true;
    } else {
        deviceIdx = uint.MaxValue;
        btnMeasSinTpx3.Enabled = false;
    }
}
}

```

```

private void btnMeasSinTpx3_Click(object sender, EventArgs e) { //
=====
    btnMeasSinTpx3.Enabled = false;

    int opm;
    if (rbOPM0.Checked) opm = 0; else opm = 2;
    int rc = pxcSetTimepix3Mode(deviceIdx, opm);
    msgToList("pxcSetTimepix3Mode", rc);

    double[] frameToaITot = new double[devPixels];
    UInt16[] frameTotEvent = new UInt16[devPixels];
    UInt32 size = devPixels;
    double t = Convert.ToDouble(txtMeasTime.Text);

    rc = pxcMeasureSingleFrameTpx3(deviceIdx, t, frameToaITot, frameTotEvent, ref
size, 0);
    msgToList("pxcMeasureSingleFrameTpx3", rc);

    btnMeasSinTpx3.Enabled = true;
}
}

```

```

        int[] imgData1 = new int[devPixels];
        int[] imgData2 = new int[devPixels];
        colorizeData(frameToaITot, imgData1);
        String colInfo1 = $"min:{colorizeMin:F}, max:{colorizeMax:F}, hit
pixels:{colorizePixCount:D}";
        colorizeData(frameTotEvent, imgData2);
        String colInfo2 = $"min:{colorizeMin:F}, max:{colorizeMax:F}, hit
pixels:{colorizePixCount:D}";
        pbFrame1.SizeMode = PictureBoxSizeMode.StretchImage;
        pbFrame2.SizeMode = PictureBoxSizeMode.StretchImage;

        IntPtr imgDataPtr1 = Marshal.AllocHGlobal(Marshal.SizeOf(typeof(int)) *
(int)devPixels);
        Marshal.Copy(imgData1, 0, imgDataPtr1, (int)devPixels);
        IntPtr imgDataPtr2 = Marshal.AllocHGlobal(Marshal.SizeOf(typeof(int)) *
(int)devPixels);
        Marshal.Copy(imgData1, 0, imgDataPtr2, (int)devPixels);

        Bitmap bmp1 = new Bitmap((int)devXsize, (int)devYsize, (int)devXsize * 4,
System.Drawing.Imaging.PixelFormat.Format32bppRgb, imgDataPtr1);
        Bitmap bmp2 = new Bitmap((int)devXsize, (int)devYsize, (int)devXsize * 4,
System.Drawing.Imaging.PixelFormat.Format32bppRgb, imgDataPtr2);

        pbFrame1.Image = bmp1;
        pbFrame2.Image = bmp2;

        labImgInfo.Text = $"Frame1 - ToA/IToT:\n{colInfo1:S}\n\nFrame2 -
ToT/Event:\n{colInfo2:S}";
    }
}
}

```

Related

- [Binary core API](#)
- [Pixet SDK](#)
- [Simple C# commandline example](#)