

# PIXet Basic

Online version: [https://wiki.advacam.cz/wiki/PIXet\\_Basic](https://wiki.advacam.cz/wiki/PIXet_Basic)

# Contents


<b>Introduction</b> .....	3
<b>Pixet Basic Interface Description</b> .....	3
<b>Title bar</b> .....	4
<b>Menu bar</b> .....	4
File .....	4
View .....	4
Settings .....	4
Scripting .....	5
<b>Quick access toolbar</b> .....	5
<b>Device list column</b> .....	6
<b>Image panel</b> .....	6
<b>Right panels</b> .....	6
Image properties .....	6
Measurement .....	6
Image Info .....	7
Histogram .....	7
Spectra .....	7
<b>Status bar</b> .....	7
<b>Getting started with the PIXet Basic</b> .....	8
<b>Collecting Data</b> .....	8
<b>Exporting data</b> .....	9
<b>Related</b> .....	10



## Introduction

The PIXet Basic is minimalist option of the [PIXet Pro](#), intended to be used only with the [MiniPIX-EDU](#) device. The PIXet Pro switches to the Basic mode when the [MiniPIX-EDU](#) is connected.

## Pixet Basic Interface Description



- 1: Title bar
- 2: Menu bar
- 3: Quick access toolbar
- 4: Device list column
- 5: Image panel
- 6: Image properties
- 7: Measurement
- 8: Image info
- 9: Spectra

## 10: Status bar

### Title bar

"Pixet Basic", {Pixet vesion}, "Build" {build\_version}, (release date), "MiniPIX", {Chip ID}, "sn" {serial\_number}

### Menu bar

#### File

---

#### Open frame...

Opens a dialog that allows to open previously recorded frame files using a \*.dsc file

#### Save data...

Opens a dialog that allows to save current frame in various file formats

Relevant file formats: .txt, .clog, .png, .tif, .tiff

#### Load factory config

Looks into (Pixet\_folder)\factory\ and searches for MiniPIX-{chip\_ID}.xml and configures the device by it.

#### View

---

#### Mirror Image

Toggle option, allows to display the detector window mirrored along Y axis

#### Rotate image

allows to rotate displayed data in detector window by 90°, 180° or 270° clockwise

#### Show grid

Toggle option, allows to display an overlay in detector window showing the pixel grid

#### Ensure aspect ratio

toggle option, keeps aspect ratio constant when resizing the Pixet window

#### Settings

---

Change password for Pixet Pro




internal development option

## Scripting

(Menu Scripting, or click on the icon with two papers at the toolbar.)

Python scripting opens Python scripting plugin, that enables to control the device using Python

template\_script.py runs an example script. If unaltered, performs 3x 1s acquisition and saves the frames into Pixet parent directory as .png files



Python plugin (Pixet Pro below it - Basic has fewer icons)

## Quick access toolbar



Quick access toolbar, left to right:

- Open frame...
- Save data...
- Show grid
- Rotate clockwise
- Color scheme
- Under warning
- Over warning
- Python scripting
- (long gap)

## Switch to PRO button

Open frame... - see File --> Open frame...

Save data... - see File --> Save data...

Show grid - see View --> Show grid

Rotate clockwise - see View --> Rotate image

### Color scheme

Allows selection of different colour representation for the measured data

### Under warning

Toggle option, highlights pixels below min value in green

### Over warning

Toggle option, highlights pixels above max value in red

### Python scripting

see scripting --> Python scripting

### Switch to PRO button

developer option, allows to switch to PRO version

## Device list column

displays all currently connected devices (their type and chip ID)

## Image panel

shows results of current measurement. Color scheme below the window can be changed using the color scheme icon. Slider at the bottom of the screen allows to browse different frames (if there are any).

## Right panels

### Image properties

---

allows to set min and max values for view of the current measurement

### Measurement



---

Mode	Imaging - measures number of hits in each pixel during selected time period Tracking - measures energy of hits in each pixel during selected time period
Frames	set how many frames will be measured
Exposure [s]	set the exposure duration of each frame
Sum	if checked, sums all measured frames in the detector window
Play	starts the measurement. Automatically repeats when finished
Image	takes a single image with selected exposure time

## Image Info

---

[X,Y]	gives X,Y coordinates of selected pixel or pixel under the mouse pointer
Count	gives number of hits in selected pixel or pixel under the mouse pointer
Min	min number of hits in image or selected area
Max	max number of hits in image or selected area
Pixel count	number of hit pixels in image or selected area
Total	total number of hits in image or selected area
Mean	mean number of hits in image or selected area

## Histogram

(only if imaging mode is selected)

## Spectra

(only if tracking mode is selected) Energy [keV]

Size [pix]

Timeline

## Status bar


Messages	shows system messages. Broom icon = clear log. Toggle option to automatically open when new message
T	duration of current experiment
A	current/total no. of measurements. Does not work when repeating or when taking single images
Measurement progress bar	Shows progress of current measurement. Does not work when repeating or when taking single images

## Getting started with the PIXet Basic


Quick intro to using the Pixet Basic program. For the Experiments cookbook, see: [EDU Kit experiments cookbook](#)

### Collecting Data

1. Connect the MiniPIX EDU camera, wait for 10 seconds and launch the software Pixet Basic.




2. The settings supposed to be modified while performing the experiments.  
See right image: The two arrows show the only place to modify the settings.





3. The collected visual data will be displayed on the black screen and the additional data (counts, histograms, etc.) is displayed in the two sections- Image Info and Spectra.

See right image: The data of the *experiment will be displayed in the two sections*



4. The image info shows energy in the “Info” tab and counts of individual particles in the “Tracks” tab.

5. The spectra section shows the energy, size, and time of individual particles in respective tabs.

Note: To enlarge the spectra results, click on the rectangle at the top right corner of the spectra window.

## Exporting data

1. Select the energy graph in spectra and right-click on the histogram in the spectra section and click on Save to File.




Figure 4. The histogram of collected data is being exported as ASCII

2. Enter the file name and select ASCII Vertical in Save as type drop-down menu. Click on Save.




Figure 5. The type of file in the save window should be vertical to make sure that we get the whole histogram

3. Next open Excel. Click on Open > Browse and browse to the saved file. Select All Files from the drop-down menu next to the File name field, select the file and click open.




Figure 6. The Open window of excel showing the "vtxt" file

4. Select Delimited in the Text Import Wizard window and click on Finish.




Figure 7. Excel delimiting the data to sort it in different columns

5. The first two columns of sheet are the list of the energy of alpha particles, the next two columns are for electrons, the next two are muons and the last two are dots.

6. Now select the first two columns and click on the Insert menu. In the Charts section, click on the Scatter chart and select Scatter with Straight Lines to get the histogram.

## Related

- [MiniPIX-EDU](#)
- [EDU Kit experiments cookbook](#)



- [The PIXet Pro program](#)

